2017-01
801515
1.0
M

 TRACOE medical GmbH
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200
info@tracoe.com, www.tracoe.com

TRACOE®
twist

CE
0 1 2 3

TRACOE® *twist* Tracheostomy Tubes




TRACOE® *twist* Tracheostomiekänäülen

**REF 301, REF 302, REF 303, REF 304
REF 305, REF 306, REF 888-306, REF 309**

EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
DA	Brugsanvisning
FI	Käyttöohje
NB	Bruksanvisning
HU	Használati utasítás
SV	Bruksanvisning
NL	Gebruiksaanwijzing
HR	Upute za uporabu
SL	Navodila za uporabo
CS	Návod k použití
RO	Instrucțiuni
TR	Kullanma talimatı
RU	инструкция по применению
EL	Οδηγίες χρήσης
KO	사용 설명서
ZH	用于新生儿
AR	دليل استعمال

TRACOE®

 TRACOE medical GmbH
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200
info@tracoe.com, www.tracoe.com

EN / Warning. Please read Instructions for Use • **DE** / Achtung. Gebrauchsanweisung beachten • **FR** / Attention. Respecter le mode d'emploi • **IT** / Attenzione. Consultare le istruzioni per l'uso • **ES** / Atención. Lea atentamente las instrucciones de uso • **PT** / Atenção. Leia as instruções de utilização • **DA** / Vigtigt. Følg brugsanvisningen • **FI** / Huomio. Noudata käyttöohjetta • **NB** / OBS. Følg bruksanvisningen • **HU** / Figyelem. Olvassa el a használati útmutatót • **SV** / OBS. Följ bruksanvisningen • **NL** / Waarschuwing. Voor gebruik de gebruiksaanwijzing lezen • **HR** / Upozorenje. Pridržavajte se uputa za uporabu • **SL** / Pozor. Upoštevajte navodila za uporabo • **CS** / Pozor. Respektujte návod k použití • **RO** / Atenție. Respectați instrucțiunile de utilizare • **TR** / Dikkat. Kullanma talimatına dikkat ediniz • **RU** / Внимание! Соблюдать инструкцию по применению • **EL** / Προσοχή. Τηρήστε τις οδηγίες χρήσης • **KO** / 경고. 사용 설명서를 읽으십시오. • **ZH** / 注意。请阅读使用说明 • **AR** / تنبيه: يجب مراعاة دليل الاستعمال

EN / For single-patient use only • **DE** / Nur für einen Patienten • **FR** / Pour un seul patient • **IT** / Prodotto destinato ad un solo paziente • **ES** / Para uso por un solo paciente • **PT** / Só para um paciente • **DA** / Kun til én patient • **FI** / Vain yhden potilaan käyttöön • **NB** / Kun for én pasient • **HU** / Csak egy betegnél alkalmazható • **SV** / Endast för en patient • **NL** / Voor gebruik bij slechts één patiënt • **HR** / Samo za jednog pacijenta • **SL** / Samo za enega bolnika • **CS** / Pouze pro jednoho pacienta • **RO** / Produs destinat unui singur pacient • **TR** / Sadece tek bir hasta için • **RU** / Только для одного пациента • **EL** / Μόνο για έναν ασθενή • **KO** / 단일 환자 전용 • **ZH** / 仅限于一位患者使用。 • **AR** / تستعمل لمريض واحد فقط

REF **EN** / Item number • **DE** / Artikelnummer • **FR** / Numéro d'article • **IT** / Codice prodotto • **ES** / Número de artículo • **PT** / Número do artigo • **DA** / Artikelnummer • **FI** / Tuotenumero • **NB** / Artikkelnnummer • **HU** / Cikkszám • **SV** / Artikelnummer • **NL** / Artikelnummer • **HR** / Broj artikla • **SL** / Številka artikla • **CS** / Číslo výrobku • **RO** / Număr articol • **TR** / Ürün numarası • **RU** / Номер артикула • **EL** / Αριθμός είδους • **KO** / 품목 번호 • **ZH** / 产品编号 • **AR** / رقم السلعة

LOT **EN** / Batch code • **DE** / Chargencode • **FR** / Numéro de lot • **IT** / Numero di lotto • **ES** / Código de lote • **PT** / Código do lote • **DA** / Batchkode • **FI** / Eräkoodi • **NB** / Batchkode • **HU** / Tételkód • **SV** / Satskod • **NL** / Lotcode • **HR** / Oznaka serije • **SL** / Koda šarže • **CS** / Kód šarže • **RO** / Număr lot de producție • **TR** / Parti kodu • **RU** / Код партии • **EL** / Κωδικός παρτίδας • **KO** / 배치 코드 • **ZH** / 批次编号 • **AR** / رمز الشحن

ONLY

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Instructions for Use for TRACOE® *twist* Tracheostomy Tubes

NOTE:

- Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own and your patients' safety, please observe the following safety information.

- The illustrations to which the text refers can be found on the (foldout) illustrated pages at the beginning of these instructions. The numbers and letters in brackets refer to the respective illustrations and product components of the tracheostomy tube. The symbols and icons used are explained on pages 1 to 5.

1. General Information

CAUTION: The product may only be used by doctors and other **trained** staff responsible for the patients' care! The attending doctor must determine the correct size and type.

Intended Purpose: This medical device is a tracheostomy tube.

Functional Description: The tracheostomy tube ensures that respired air passes through the tracheostoma.

- The inflated cuff separates the upper from the lower respiratory tract, preventing the flow of air from the lungs into the mouth and nose and vice versa. This ensures that the patient can breathe **only** through the tube. The cuff is inflated by inflating it with air using the inflation line. When inflated, the cuff seals the space between the trachea and the outer cannula wall. This allows efficient ventilation with ventilation equipment and prevents subglottic secretions from entering the lower respiratory tract.
- Fenestrated cannulas allow some of the air needed for speaking to reach the upper respiratory tract through the fenestration. Airway resistance in the upper respiratory tract is also reduced.
- See below for a more detailed functional description.

Maximum Useful Life of the Tracheostomy Tube:

- With cuff – 29 days,
 - Without cuff – 5 weeks,
- from first use (see section 5 "General Precautions"). This maximum period of use also includes any periods when the tracheostomy tubes or inner cannulas are not in use.

Product for Single-Patient Use: Multiple use by one and the same patient during the period of use is permitted.

Product Info Card: This product is supplied with a product info card with two detachable labels. The labels show product details. Keep the info card separately to facilitate reordering. The detachable label can be affixed to the patient's medical records, for example.

2. General Description

The product consists of an outer cannula (1) with a cuff (2) or without a cuff, and an exchangeable inner cannula with a 15-mm connector (3) or grooved low profile connector .

(1) and (2) are connected and disconnected with a twist lock (C). The cuff is inflated or deflated using the inflation line (2a).

The outer cannula is made of radiopaque polyurethane and is attached to a neck flange (5) that can swivel on 2 axes (cardanic suspension).

The perforated obturator (6), in sizes 06 – 10, allows the use of a Seldinger wire ($\varnothing 1.27 \pm 0.04$ mm) for cannulation.

3. Indications

The products are suitable for patients where access to the respiratory tract is required via a tracheostoma.

- Tubes without a cuff are particularly suited for patients who breathe spontaneously but need to use a tube.
- Tubes with a cuff are suitable for patients where access to the respiratory tract is required via a tracheostoma with the trachea sealed.
- Fenestrated tubes facilitate speech in patients with an intact larynx.
- Tubes with a suction line are for use in patients where suction of subglottic secretion is indicated.
- REF 309 tubes with an air supply for speaking allow speech in patients with an intact larynx who are awake and mechanically ventilated.
- The REF 305 tube is especially suitable for patients in whom a tube is indicated after laryngectomy.

4. Contraindications

4.1 Absolute contraindications

- Do not use the occlusion cap/speaking valve in laryngectomised patients (patients without a larynx) – risk of suffocation!

4.2 Relative contraindications

Relative contraindications are cases where the risks and benefits of the procedure must be carefully evaluated:

- Paediatric use.
- Unusually deep trachea, e.g. in obese patients, where a longer tube may be required.

5. General Precautions

- You are strongly advised to keep a ready-to-use spare tube and several spare inner cannulas at the patient's bedside. Store in a clean and dry state.
- Check the tube for integrity and correct function prior to use/insertion. Check, for example, for the absence of obstruction, cuff seal, correct and stable fit of the inner cannula inside the outer cannula, absence of kinks, stable connection between the tube and neck flange, etc. The cuff material must not be brittle. If the product is damaged, replace with a new product.
- Do not use force on the tracheostomy tube, as this may damage or break it. If the connections at the 15 mm connector are tightly fixed, always use a disconnect wedge approved for tracheostomy tubes.
- Keep the 15 mm connector **(3)** clean and dry.
- When changing the inner cannula, always ensure that the inflation line **(2a)** of the cuff is not positioned between the inner and outer cannulas, as it may get trapped and damaged.
- When repositioning the patient in bed, ensure that the patient does not lie on the pilot balloon **(2a)**, as this would greatly increase the cuff pressure and may damage the trachea.
- During mechanical ventilation and frequent changes of the patient's position or manipulation of the tube, the inner cannula may become separated from the outer cannula.
- With the use of high ventilation pressures, leakage between the inner and outer cannulas may occur in isolated cases.
- The use of tracheostomy tubes can lead to pressure points, necrosis of the skin on the neck and skin irritation (e.g. moisture). To avoid this, we recommend placing a pad underneath the neck flange.
- To avoid damage to the material, the cuff should not be in contact with lidocaine-containing aerosols or any ointments.
- All components of the cuff inflation system must be kink-free and not taut when testing the cuff pressure. The hand-held manometer may otherwise give a false pressure reading.
- Thin cuffs are to some extent permeable to water vapour. Therefore, condensate may accumulate in the cuff. Where quantities are small, this is of no significance. However, if larger

quantities of condensate are inadvertently aspirated into the inflation line, it is no longer possible to measure and adjust the cuff pressure correctly and the tube needs to be replaced. Before removing the tube, remove as much air or water from the cuff as possible using a syringe.

- When using together with other medical devices, the respective instructions for use must be observed. If in doubt, contact the manufacturer.

6. Warnings

- Do not use if the sterile packaging is not intact.
- Never use fenestrated inner cannulas for ventilation.
- Use the occlusion cap only with fenestrated tubes (outer + inner cannula) and an unblocked cuff.
- Do not use a speaking valve during sleep.
- When using laser or electrosurgical treatment equipment, ensure that the tracheostomy tube is at a safe distance. There is a risk of fire and toxic gases and the tube may get damaged.
- Select the correct tube to ensure that the fenestration (if applicable) is positioned at sufficient distance from the stoma canal. Failure to do so exposes ventilated patients to the risk of emphysema. Granulation tissue may also form or airway resistance may increase if speaking valves or occlusion caps are used.
- If a fenestrated tube is used, there is a risk of increased formation of granulation tissue.
- The cuff pressure can change, for example due to altitude changes (e.g. in an aircraft), or if laughing gas is used as an anaesthetic or during the connection and disconnection of a hand-held manometer.
- Excessive cuff pressure poses a risk of permanent damage to the trachea.
- Excessive cuff pressure may give rise to cuff hernias.
- Insufficient cuff pressure poses a risk of aspiration.
- Ensure that the correct Luer connectors are used (e.g. **2a + 9**).
- During insertion and removal of the tube, irritation, need to cough or bleeding may occur.

7. Adverse Reactions

Pressure points, necrosis, skin irritation, granulation tissue, need to cough, difficulty swallowing, and bleeding.

8. Inserting the Tube

8.1 Preparing the Tube

1. Check the pack content for completeness (**D**).
2. For tubes with a cuff, check the cuff (**2**) for leakage by inflating it. For this purpose, inflate the cuff with a hand-held manometer to a pressure of 50 cmH₂O (\approx 36.78 mmHg) and wait for 1 minute to see whether the cuff becomes deflated. If the cuff is leak-tight, remove all air from the cuff with a syringe. Push the cuff upwards in the direction of the neck flange (**5**) to facilitate insertion.

8.2 Preparing the Patient

- Ensure that the patient is well pre-oxygenated before cannulation or recannulation.
- Slightly overextend the patient's neck, if possible, to facilitate insertion.
- Should complications arise during cannulation or recannulation, take appropriate safety measures (e.g. tracheal spreader, tube with a smaller diameter) to allow the doctor to carry out short-term ventilation using translaryngeal intubation or a larynx mask.

8.3 Inserting the Tube

Carry out the following steps:

1. Before inserting the tube, remove the inner cannula (**3**) and fit the obturator (**6**) into the outer cannula (**1**). When inserting the tube, hold the tube by the flange (**5**) and press the obturator firmly into the upper end of the outer cannula with the thumb. To aid insertion, a thin layer of water-soluble lubricant may be applied to the protruding part of the obturator at the patient end and to this area of the tube including the cuff.
2. Once the tube has been inserted into the patient's trachea, remove the obturator immediately.
3. Then insert the inner cannula. To secure the inner cannula, hold the neck flange with the fingertips and twist the connector of the inner cannula until it locks into place, i.e. until the blue arrows and/or markings are aligned with each other (**C**).
4. Attach the neck strap (**12**) to the neck flange (**5**) to keep the tube in place.
5. Check the position of the tube as well as its correct function (see chapter 6 "Warnings"). If required, carefully suction the trachea and tube to ensure that the airways are unobstructed.
6. For ventilation, securely connect the standardised 15 mm connector (**3**) of the inner cannula to a ventilation machine.

For cuffed tubes:

7. Inflate the cuff using the pilot balloon (**2a**). Now check the cuff pressure to ensure that the cuff has not been damaged

during insertion. Adjust the cuff pressure to the individual ventilation therapy and check at regular intervals. Typically, the pressure should be between 20 cmH₂O (\approx 15 mmHg) and 30 cmH₂O (\approx 22 mmHg). Attaching a hand-held manometer causes a fall in cuff pressure. This effect is more marked with small sizes.

8.4 Removing the Tube (see chapter 8.2)

Prepare to remove the tube as follows:

- Bend the patient's head slightly backwards and
- completely unblock the cuff of cuffed tubes (see chapter 8.5).

8.5 Unblocking the Cuff

Before unblocking, ensure that as little secretion as possible can enter the bronchi.

While the cuff is being unblocked, suction secretions with a suction catheter passed through the tube. If tubes with subglottic suction have been used, also aspirate the subglottic area before unblocking the cuff, see chapter 9.5.

To unblock the cuff, connect a syringe to the pilot balloon (2a). Completely remove all air/condensate if applicable. (See chapter 5. "General Precautions").

If it is not possible to unblock the cuff, apply the following measures:

- Remove the inner cannula. Make sure that the inflation line (2a) is not kinked and unblock the cuff again.
- If it is still not possible to unblock the cuff: carefully cut through the inflation line between the 15 mm connector (3) and the neck flange (5).

9. Handling

9.1 Changing the Inner Cannulas

If viscous secretion collects in the inner cannula and cannot be suctioned, thus impeding the airflow, replace the inner cannula with a new or cleaned inner cannula.

To remove the inner cannula, turn the 15 mm connector of the inner cannula anticlockwise (C).

Once a new inner cannula has been inserted into the outer cannula, turn the 15 mm connector of the inner cannula clockwise until it locks into place, i.e. until the blue arrows and/or embossed markings are aligned with each other (C).

Inner cannulas with grooved low profile connectors function in the same way.

When changing the inner cannula, always ensure that the inflation line (**2a**) of the cuff is not positioned between the inner and outer cannulas, as it may get trapped and damaged.

9.2 Fenestrated Tracheostomy Tubes

Unfenestrated inner cannulas have white 15 mm connectors and are used, for example, during ventilation.

To allow the patient to speak, insert a **fenestrated inner cannula** (blue 15 mm connector or blue grooved low profile connector) **in the fenestrated outer cannula**. Once the cuff has been unblocked, a speaking valve can be placed on the 15 mm connector of this inner cannula. Follow the instructions for use of the respective speaking valve.

The inner cannula REF 401 (not supplied) has an integrated speaking valve. The valve allows air to enter the trachea through the tube on inspiration; on expiration, the valve closes and the air flows alongside the tube, as well as through the fenestration, through the larynx into the upper respiratory tract, which allows vocalisation.

If there is an urge to cough or possible breathing difficulties, the REF 401 valve leaflets can be folded back and the patient can breathe unimpeded.

CAUTION: Speaking valves must only be used in patients who are awake and can breathe spontaneously. When a speaking valve is introduced, the patient must be instructed by trained personnel and carefully monitored while becoming accustomed to it; make sure that the patient can breathe sufficiently.

Unblock the cuff before using a speaking valve. Do not use a speaking valve during sleep.

CAUTION: Patients with laryngeal stenosis, vocal cord paralysis, severe tracheal stenosis, respiratory tract obstructions, respiratory tract infections or copious mucous secretions must not use a speaking valve.

9.3 Occlusion Cap of the Fenestrated Tube for Weaning

The fenestrated inner cannula can be closed by attaching the enclosed occlusion cap (**10**) to the 15 mm connector.

CAUTION: When preparing the tube, ensure that the patient's upper respiratory tract is unobstructed. The upper respiratory tract can be cleared by coughing out or suctioning the secretions. The

cuff needs to be unblocked. When closing the tube, the patient's breathing and vital signs must be monitored by a nursing professional. If any signs of breathing difficulties are detected, the occlusion cap must be removed immediately.

9.4 Keeping the Fenestration Open

If the tube is used over a prolonged period, it must be ensured at regular intervals that the fenestration is not blocked by secretions, incrustation or in-growing tissue. If required, replace the tube.

9.5 Application of REF 306, REF 888-306 with Subglottic Suction Line

These tubes are fitted with a flat suction line (9) that is fixed to the outer curve of the outer cannula and ends with two openings immediately above the cuff (2). The suction line has a female Luer connector, through which any secretion can be removed using a syringe. As an alternative, a special suction device with a vacuum controller connected to the enclosed connectors (11) can be used. After suctioning, ensure that the Luer connector is closed.

CAUTION:

- Ensure during suctioning that no excessive vacuum is used over a prolonged period (- maximum of 200 mbar).
- To prevent the subglottic area from drying out, we recommend the use of intermittent suction.
- The suction line may get blocked due to accumulated secretions or adhesions to the tracheal mucosa. If it is intended to **rinse the suction line** (e.g. with a small amount of air or physiological saline solution), ensure beforehand that the cuff is sufficiently blocked (risk of aspiration).
- The suction line may cause pressure points in the area of the tracheostoma or may cause the tracheostoma to become irregular. In such cases, the doctor needs to decide whether this type of tube can still be used.

9.6 Use of REF 309 with air supply for speaking

This tube is equipped with an air supply line that ends in the middle of the outer curve of the outer cannula, and has four openings. The air supply line is fitted with a T connector. At the conical opening, the line can be connected to the air or oxygen supply on the intensive care ward. Opening the side outlet allows the air or oxygen supplied to exit. With deliberate closure of the side opening, the air or oxygen flows through the trachea into the larynx, which allows vocalisation.

CAUTION: The air or oxygen supplied for speaking must always be humidified: the quantity supplied must not exceed 8 – 12 l/min.

10. Cleaning, Disinfection and Storage

10.1 Cleaning

The tube (1 + 3) needs to be cleaned before reuse. Lukewarm potable water can be used for cleaning. TRACOE medical recommends the use of the cleaning products offered by TRACOE for cleaning. Once the tubes have been cleaned, rinse with potable water and air-dry. When cleaning the outer cannulas, take care not to damage the cuff.

CAUTION:

- Never use aggressive domestic cleaning agents, cleaning agents for dentures or solvents with a high alcohol content, for example, as these may have a lasting adverse effect on the function.
- Do not heat the tubes to temperatures above 65 °C, as the safety of the products can then no longer be guaranteed.

10.2 Disinfection

Should disinfection appear necessary TRACOE medical can provide a list of tested disinfectants on request. This list is also available at www.tracoe.com. The use of unsuitable disinfectants can result in damage to the product.

Once the tubes have been disinfected, rinse with sterile water and air-dry.

10.3 Storage

Store the clean tube in a dry location.

11. Product Modifications

Modifications to the TRACOE products may only be carried out by employees or agents of TRACOE medical GmbH.

12. Disposal

These products may only be disposed of in accordance with the valid national regulations governing waste materials.

13. Returns

Returned used products will only be accepted subject to prior consultation and if a completed decontamination certificate is enclosed. This form is available either directly from TRACOE medical or via the website www.tracoe.com

14. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be effected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on the website www.tracoe.com.

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