TRACOE medical GmbH

Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany

Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200

info@tracoe.com, www.tracoe.com

TRACOE®
vario

Instructions for Use for

TRACOE® *vario* Tracheostomy Tubes

REF 450, REF 451, REF 455, REF 460, REF 461,
REF 464, REF 470, REF 471

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TRACOE® *vario* Tracheostomiekanülen

REF 450, REF 451, REF 455, REF 460, REF 461,
REF 464, REF 470, REF 471

2		Symbols, Symbole
6	EN	English
16	DE	Deutsch
27	FR	Français
38	IT	Italiano
49	ES	Español
60	PT	Português
71	DA	Dansk
81	FI	Suomi
91	NB	Norsk
101	HU	Magyar
111	SV	Svenska
121	NL	Nederlands
132	HR	Hrvatska
141	EL	Ελληνικά
153	KO	한국(어)
161	ZH	中文
169	AR	عربي
177		Techn. Drawing, techn. Zeichnung Sizes, Größen

EN / Warning. Please read Instructions for Use · **DE** / Achtung. Gebrauchsanweisung beachten · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter atenção às instruções de utilização · **DA** / Vigtigt: Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NB** / OBS. Følg bruksanvisningen · **HU** / Figyelem. Olvassa el a használati útmutatót · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **HR** / Upozorenje. Pridržavajte se uputa za uporabu · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오. · **ZH** / 注意。注意使用说明书 · **AR** / تنبيه. يجب مراعاة دليل الاستعمال

EN / For single patient use only · **DE** / Nur für einen Patienten · **FR** / Pour un seul patient · **IT** / Prodotto destinato ad un solo paziente · **ES** / Sólo para un paciente · **PT** / Só para um paciente · **DA** / Må kun bruges til én patient · **FI** / Vain yhden potilaan käyttöön · **NB** / Kun for én pasient · **HU** / Csak egy betegnél alkalmazható · **SV** / Endast för en patient · **NL** / Slechts voor één patiënt bestemd · **HR** / Samo za jednog pacijenta · **EL** / Μόνο για έναν ασθενή · **KO** / 오직 일인 환자 용입니다 · **ZH** / 只供一位病人使用 · **AR** / تستعمل لمريض واحد فقط

REF **EN** / Item number · **DE** / Artikelnummer · **FR** / Référence · **IT** / Codice prodotto · **ES** / Número de artículo · **PT** / Número do artigo · **DA** / Artikelnummer · **FI** / Tuotenumero · **NB** / Artikkelnummer · **HU** / Cikkszám · **SV** / Katalognummer · **NL** / Artikelnummer · **HR** / Broj artikla · **EL** / Αριθμός είδους · **KO** / 품목 번호 · **ZH** / 产品编号 · **AR** / رقم السلعة

LOT **EN** / Batch code · **DE** / Chargencode · **FR** / Code de lot · **IT** / Numero di lotto · **ES** / Código de lote · **PT** / Código do lote · **DA** / Batchkode · **FI** / Eräkkoodi · **NB** / Chargekode · **HU** / Tételkód · **SV** / Satsnummer · **NL** / Lotcode · **HR** / Oznaka serije · **EL** / Κωδικός παρτίδας · **KO** / 생산 단위 · **ZH** / 批次編碼 · **AR** / رمز الشحن

EN / Date of manufacture · **DE** / Herstellungsdatum · **FR** / Date de fabrication · **IT** / Data di produzione · **ES** / Fecha de fabricación · **PT** / Data de fabricação · **DA** / Produktionsdato · **FI** / Valmistuspäivä · **NB** / Produksjonsdato · **HU** / Gyártási dátum · **SV** / Tillverkningsdatum · **NL** / Fabricagedatum · **HR** / Datum proizvodnje · **EL** / Ημερομηνία κατασκευής · **KO** / 제조일 · **ZH** / 生产日期 · تاريخ التصنيع / **AR**

EN / Manufacturer · **DE** / Hersteller · **FR** / Fabricant · **IT** / Produttore · **ES** / Fabricante · **PT** / Fabricante · **DA** / Producent · **FI** / Valmistaja · **NB** / Produsent · **HU** / Gyártó · **SV** / Tillverkare · **NL** / Fabrikant · **HR** / Proizvođač · **EL** / Κατασκευαστής · **KO** / 제조사 · **ZH** / 生产商 · الشركة المصنعة / **AR**

EN / Use by date · **DE** / Verwendbar bis · **FR** / Utilisable jusqu'au · **IT** / Utilizzare entro il · **ES** / Fecha de caducidad · **PT** / Utilizável até · **DA** / Anvendes inden · **FI** / Käytettävä ennen · **NB** / Kan brukes til · **HU** / Felhasználhatóság dátuma · **SV** / Används före-utgångsdatum · **NL** / Bruikbaar tot · **HR** / Upotrijebiti do · **EL** / Χρήση έως · **KO** / 유통기한 · **ZH** / 保质期至 · قابل للاستخدام حتى / **AR**

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Instructions for Use for TRACOE® *vario* Tracheostomy Tubes

These instructions only apply to the following TRACOE® *vario* tracheostomy tubes:

- REF 450** Tracheostomy tube with low-pressure cuff, spiral reinforcement and adjustable flange
- REF 451** Extra-long tracheostomy tube with low-pressure cuff, spiral reinforcement and adjustable flange
- REF 455** Tracheostomy tube with spiral reinforcement and adjustable flange
- REF 460** Tracheostomy tube with low-pressure cuff and adjustable flange
- REF 461** Extra-long tracheostomy tube with low-pressure cuff and adjustable flange
- REF 464** Tracheostomy tube, fenestrated with adjustable neck flange
- REF 470** Tracheostomy tube with low-pressure cuff, adjustable flange and suction line
- REF 471** Extra-long tracheostomy tube with low-pressure cuff, adjustable flange and suction line

1. General Information

CAUTION: The following instructions must be read carefully before using the products. These instructions exclusively refer to

the products listed above. Individuals using the products must be familiar with these instructions.

Sterilisation: Contents are sterile unless the seal is broken or the package is opened or damaged. The products must not be sterilised again after use.

Single-Patient Product: TRACOE® *vario* tracheostomy tubes are medical products intended for use with a single patient only. Repeated use with the same patient is permitted.

Useful Life: The manufacturer recommends that TRACOE® *vario* tracheostomy tubes with and without cuff should not be used for more than 29 days. Even if tubes are slightly damaged, they must absolutely be replaced. Frequent replacement of tracheostomy tubes is recommended and should be supervised by the responsible doctor.

Composition/Disposal: TRACOE® *vario* tracheostomy tubes are made of materials allowed internationally for medical use. They do not contain latex. Disposal of these products should comply with the relevant national waste disposal regulations.

Pack Contents: Each pack contains one tube, one obturator to assist insertion, a neck strap, one instruction manual and a tube information card. The REF 470 and REF 471 tube also comes with connectors for the suction tube of a suction pump.

Sizes: TRACOE® *vario* tubes REF 450, REF 455 and REF 460 are available in 6 sizes (sizes 06 to 11). The REF 451, REF 461, REF 470 and REF 471 tubes are available in sizes 07 to 10. REF 464 is available in sizes 07, 08 and 09. The size refers to the internal diameter of the tracheostomy tube.

Dimensions: The dimensions are shown in table 1. Length B indicates for each tube size the distance over which the adjustable flange at the machine end of the tube can be adjusted.

Product Info Card: This product is supplied with a product info card with two detachable labels. The labels show product-specific data. The info card should be kept separately as it makes reordering easier. The detachable label can be attached, e.g., to the patient's record.

2. Description

TRACOE® *vario* tracheostomy tubes consist of a flexible tube with an obturator and a neck flange which is movable lengthwise with respect to the tube.

2.1 Tubes

The tubes are bent at $100 \pm 5^\circ$ and have a 15 mm connector at the proximal (machine) end. The tubes REF 450, REF 451 and REF 455 are reinforced by a spiral that is fully embedded in the tube wall and can therefore not be bent or squashed. The REF 450, REF 460, REF 461, REF 470 and REF 471 tubes have a low-pressure cuff at the distal (patient) end of the tube.

Tubes REF 460, REF 461, REF 470 and REF 471 do not contain any spiral reinforcement.

The fenestrated non-reinforced REF 464 cannula allows part of the air which is needed for speaking to reach the upper airways through the fenestration. In addition, the airway resistance in the upper airways is reduced. The fenestration of the outer cannula is dimensioned so that suction tubes cannot penetrate through the holes into the trachea.

REF 470 and REF 471 have a suction line integrated in the wall of the tube with an opening above the cuff and a suction tube at the proximal (machine) end.

All TRACOE® *vario* tubes have a scale behind the 15 mm connector on its cranial surface, that is, the surface towards the connector. This scale is used by nursing staff to determine the depth to which the tube is inserted into the tracheostoma. At the same time the scale indicates the position of the tube relative to the neck flange.

CAUTION: It is essential that the scale is always oriented cranially, that is, towards the chin.

Radiopaque

Tubes with spiral reinforcement (REF 450, REF 451 and REF 455) or X-ray contrast lines (REF 460, REF 461, REF 464, REF 470 and REF 471) are easy to see on the X-ray image.

MRT Imaging

TRACOE® *vario* with spiral reinforcement (REF 450, REF 451 and REF 455) are ferromagnetic and heat up in magnetic fields. Therefore, they cannot be used in MRI. This also affects examinations on more distant examination fields, because the tubes may be drawn into the magnetic field.

TRACOE® *vario* tubes without spiral reinforcement (REF 460, REF 461, REF 464, REF 470 and REF 471) can be used for the MRI. It is recommended that the pilot balloon with the valve be fixed outside the examination field using tape.

2.2 Adjustable Flange

All TRACOE® *vario* tracheostomy tubes have an adjustable flange at the proximal end behind the 15 mm connector. The flange is first positioned by means of a spring mechanism and, in the next step, fixed in place by moving the orange lever. Before this locking lever is used, the flange can be moved along the tube by pressing the spring element, without changing the position of the tube in the trachea. Releasing the spring leaves the flange in its new position. In this way, the flange can be moved easily by pressing it with two fingers. The ability to move the flange is of great importance for the individual adjustment of the tube length at the neck, for example in the event of postoperative oedema or for patients with a thick neck.

Once the flange is in the desired position, the orange locking lever on the right-hand side of the flange must be pushed upwards  to operate the locking device and fix the flange in place. If it is pushed downwards  the locking function is released. When the locking device is released, the flange can be moved along the tube by pressing the spring element as described above and be brought into position.

The flange includes two flexible wings bent at an angle of approx. 45° to the patient with eyelets at the end for fastening the neck strap. The wings can be turned together around the flange, or rotated independently within a range of about 60° to 180°. This permits universal adjustment to any particular neck anatomy.

3. Indications

The products are suitable for all patients where access to the airways via a tracheostoma is required. The fenestrated tube REF 464 facilitates the weaning of the patient from the tube and, in combination with a speaking valve, promotes speaking in patients with an intact larynx. All speaking valves allow the air to flow through the tube into the trachea during inhalation. When the patient exhales, the valves close and the air flows through the *larynx* into the upper airways, both past the tube and through the fenestration, thereby allowing the speech to be formed.

The TRACOE® *vario* tracheostomy tubes with cuff (REF 450, REF 451, REF 460, REF 461, REF 470 and REF 471) are suitable

for all cases where an additional sealing of the trachea is required.

The REF 470 and REF 471 tubes with low-pressure cuff and suction line are predominantly used for patients in whom suctioning of the subglottic space is indicated. The tubes with spiral reinforcement, REF 450, REF 451 and REF 455, are predominantly used in the event of narrowing processes in the trachea. The REF 451 extra-long tubes reinforced with a spiral or the REF 461 and REF 471 extra-long tubes without spiral reinforcement are indicated wherever the standard tube length is insufficient, for example in the case of a thick or swollen neck, or when it is necessary (for example with stenosis) to move the low-pressure cuff towards caudal.

The TRACOE® *vario* tracheostomy tubes without spiral reinforcement, with cuff (REF 460 and REF 461), are mostly used postoperatively with spontaneously breathing patients, or in the case of narrowing processes in the trachea, oedema in the neck or thick neck.

4. Directions for Use

4.1 Preparation of the Tube

Tube Selection

In order to guarantee optimal respiration, the largest tube that is compatible with the patient's anatomy should be selected. Determination of the size of tube to be used should be carried out by a doctor.

The article number (REF), size, internal diameter (ID), external diameter (AD) and length are all indicated on the flange of TRACOE® *vario* tracheostomy tubes so that all relevant information about a tracheostomy tube in place is easily found.

Cuff Leakage Test

The cuff is checked for leaks by test inflation. The test inflation involves inflating the cuff using a hand-held manometer, e.g. the TRACOE® cuff pressure monitor REF 720, to a pressure of 50 cm H₂O and observing whether the cuff deflates over a period of one minute. If no leakage occurs, the cuff is to be deflated completely and to be moved in the direction of the neck plate. This facilitates the subsequent insertion of the tube.

CAUTION: All parts of the cuff inflation system must be free from strain and bending during the application of the hand-held manometer, otherwise the hand-held manometer may show incorrect pressure values.

CAUTION: When a pressure measuring unit such as the TRACOE cpm is connected to the valve of a tracheostomy tube with cuff, there will always be a small amount of air flowing from the cuff into the unit's measuring chamber, which was previously depressurised. As a result, there will be a slight pressure loss in the cuff. This means that the cuff pressure before the connection of the measuring unit was always slightly higher than the pressure value that is displayed at the connected measuring unit. The pressure loss due to the filling of a measuring chamber through the cuff is normal. If the measuring unit is frequently connected without the cuff being refilled, one may get the impression that the cuff is leaking – especially for smaller tube sizes.

4.2 Inserting the Tube

The obturator must be inserted into the tube prior to insertion. When inserting the tube into the tracheostoma, the tube should be held at the flange, and the obturator pressed firmly into the 15 mm connector with the thumb. To assist insertion, a thin film of a water-soluble lubricant can be applied to the protruding part of the obturator, and the lower part of the tube including the cuff.

CAUTION: When inserting the tube, it is essential to ensure that the cuff is not punctured by instruments or sharp cartilage ridges.

Immediately after inserting the tube into the trachea, the obturator must be removed, the cuff inflated and the airway through the tube checked to ensure it is free of obstruction.

4.3 Inflating the Low Pressure Cuff

The low pressure cuff is inflated by the Luer valve located at the pilot balloon. Instead of a standard syringe we recommend the use of a hand-held manometer (e.g. REF 720). The pressure must be adjusted for the individual patient and should be between 20 cm H₂O (≈ 15 mm Hg) and 30 cm H₂O (≈ 22 mm Hg). To avoid damage to the trachea, the cuff pressure must not be higher than 30 mbar (≈ 22.5 mm Hg / ≈ 30.59 cm H₂O). The inner cuff pressure must be controlled regularly by means of a manometer to avoid damage to the trachea. Please note that, when connecting the manometer to the tube's pilot balloon, there will be a brief and slight loss of air from the tube's cuff into the measuring chamber of the hand-held manometer, which will result in a slight pressure drop inside the cuff. This is a normal physical process. If, at the maximum tolerable pressure of 30 mbar (≈ 22.5 mm Hg / ≈ 30.59 cm H₂O), a seal between the trachea and the cuff is not obtained, the cuff must be fully deflated and then re-inflated. If there

is still no seal after this second attempt, a bigger tube must be used. In case of mechanical ventilation, the correct occlusion volume is obtained when, upon a brief increase of pressure on the chest, a small amount of air passes the cuff. The same applies to spontaneous breathing if, with the proximal end of the tube occluded and during forced exhalation, a small amount of air passes between cuff and tracheal wall, audible at the open mouth.

4.4 Tube Attachment

The tube is held in place using a neck strap which passes through the two eyelets on the flange wings. These wings should be positioned according to the patient's anatomy. It is recommended that a compress is placed between the tracheostoma and the adjustable flange in order to prevent irritation of the skin underneath the flange. While attending to the tracheostoma or if oedema develops in the neck area, the flange can be adjusted as required, as described in section 2.2.

4.5 Suctioning of Secretions with Tube REF 470 resp. REF 471

If using the REF 470 resp. REF 471 tube, secretions which collect in the subglottic space can be suctioned through the opening of a channel integrated into the tube wall above the low-pressure cuff. Suctioning can be performed by attaching a syringe with a male Luer connector to the female Luer connector on the suction tube, or by using the connectors supplied to connect to the suction tube of a suction pump. If the suction tube becomes blocked, it can be cleared by blowing air through or cleaning with 2 to 3 ml of sterile saline solution.

CAUTION: When suctioning, ensure that low pressure is not applied for an extended period and/or that the subglottic space does not dry out.

4.6 Deflation of the Cuff / Removal of the Tube

Before deflating the cuff, any secretion which has collected above the cuff should be suctioned. Then the cuff should be deflated using a syringe so that the tube can be removed through the tracheostoma with minimum resistance.

4.7 Fenestrated tracheostomy tubes (REF 464)

Speaking

For speaking, a speaking valve can be attached to the 15 mm connector. The instructions for use of the respective speaking valve must be observed.

CAUTION: The fenestration must always be positioned in the trachea in a way that it is not obstructed.

CAUTION: The speaking valves may only be inserted when the patient is awake because the patient must be able to breathe spontaneously after the air has been removed from the cuff. The patients must be thoroughly instructed and monitored by trained staff in order to ensure sufficient breathing.

Weaning

For weaning, an occlusion cap for 15 mm connectors can be attached to the 15 mm connector.

CAUTION: When preparing the tube, it has to be ensured that the upper airways of the patient are clear. The upper airways are cleared by coughing up the secretion or removing it by means of suction. When closing the tube, the patient's breathing and vital signs must be monitored by the nursing staff. If the patient shows any signs of breathing difficulties, the occlusion cap must be removed immediately.

Keeping the Fenestration Open

If the tube is used over a long period of time, it has to be ensured at regular intervals that the fenestration is not blocked by secretion, incrustation or in-growing tissue. If required, the tube must be changed.

4.8 General Precautionary Measures

CAUTION: Damage Resulting from the Use of Force

Do not use force on the tracheostomy tube, because otherwise there is the risk of damage (e.g. breaking). In case of firmly attached connections that can only be removed with difficulty, a disconnecting aid that is approved for tracheostomy tubes must be used.

CAUTION: Water in Cuff

Thin cuffs have a certain degree of permeability to water vapour. For this reason, condensed water vapour may accumulate in the cuff. When quantities are small, this is of no significance. If, however, in the event of larger quantities, water is inadvertently aspirated into the inflation line, it is no longer possible to measure and adjust the cuff pressure correctly, and the tube has to be replaced. Before removing the tube, any air or water located in the cuff must be removed to the best possible extent by using a syringe.

CAUTION: Flights

Under normal conditions, we recommend that the cuff pressure is checked at least three times a day and adjusted if necessary. In particular, when being subjected to changes in air pressure, for example due to changes in altitude (e. g. in an aircraft), care must be taken to ensure that the cuff pressure is adjusted with sufficient frequency. In an aircraft or when travelling in a mountainous area, the cuff pressure may change substantially within a few minutes with serious consequences!!! In any case, the attending doctor should be consulted before taking a flight.

CAUTION:

- Never use fenestrated inner cannulas for ventilation.
- Select the tube so that the fenestration (if applicable) is positioned at sufficient distance to the stoma canal. If the above warning is not observed, there is a risk of emphysema for ventilated patients. Also, granulation tissue may form or the airway resistance may increase if speaking valves or occlusion caps are used.
- If a fenestrated tube is used, there is a risk of increased formation of granulation tissue.

5. Care, Cleaning and Disinfection of TRACOE® vario Tracheostomy Tubes

TRACOE® *vario* tracheostomy tubes are best cared for by regular and careful cleaning. We recommend the use of a complete TRACOE® tube clean set (REF 930-A or -B with cleaning powder, or REF 931-A or -B with liquid cleaning concentrate). These sets comprise a cleaning box with sieve insert, cleaning powder or liquid, swabs and a brush.

After cleaning, the tube should be rinsed with tap water or distilled water and dried. For further details see the instructions for TRACOE® tube clean sets.

CAUTION: The tubes must never be cleaned with aggressive household detergents, denture cleansers, concentrated alcohol or brushes with sharp tips as all of these may have a lasting adverse effect on how the products function. Avoid heating-up, boiling or sterilising the products with steam.

Disinfection: Since the upper respiratory tract is never free from germs – even in healthy individuals – we do not recommend the use of disinfectants.

If, however, disinfection appears necessary, TRACOE medical can provide a list of tested disinfectants on request. This list can also be accessed on the website of TRACOE medical GmbH (www.tracoe.com).

CAUTION: Under no circumstances must disinfectants be used that release chlorine or substances that contain strong alkalis or phenol derivatives.

CAUTION: The tubes should not be heated to a temperature above 50 °C, since from that temperature the safety of the products cannot be guaranteed.

CAUTION: Before re-using or re-inserting a tube, it must be checked to ensure that it is intact and in perfect working order, e.g. patency, cuff leakage, secure connection between tube and neck plate, etc.

Recommendation: We strongly recommend that you keep a spare tube in stock at all times. These tubes must be stored under clean and dry conditions.

6. Returns

Returns of used products are only accepted if they are accompanied by a filled-in decontamination certificate. You can obtain the required form either directly from TRACOE medical or from the website www.tracoe.com.

7. General Terms and Conditions

Sales, deliveries and returns of all TRACOE® products are carried out exclusively according to the General Terms and Conditions which can be obtained directly from TRACOE medical GmbH or from the website www.tracoe.com.

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