TRACOE® phon assist I Speaking Valves

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EN The warning label must be applied to the connector of the speaking valve or the inflation line of the tracheostomy tube before the speaking valve phon assist I is used.


Do **NOT** use speaking valve with **blocked cuff!**



Pls. fold
here

Do **NOT** use speaking valve with **blocked cuff!**

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Instructions for Use for TRACOE phon assist I Speaking Valves

1. General Information

TRACOE phon assist I valves are adjustable speaking valves developed for tracheostomy patients. The valves are manufactured with or without an oxygen supply port for TRACOE grid buttons and stoma buttons, TRACOE adhesive carriers and for tracheostomy tubes with a 15 mm connector.

CAUTION: the following instructions are for the use of doctors, nursing staff and the patient. The instructions refer exclusively to TRACOE phon assist I valves and must be read carefully before use.

Sterilisation: TRACOE phon assist I valves are sterilised with ethylene oxide. The contents are sterile only if the packaging and all seals are intact.

Materials: the materials used in the manufacture of TRACOE phon assist I satisfy all international requirements for medical devices; they contain no latex.

Packaging Contents: each pack REF 650-S /-SO /-T or /-TO /-TO-C contains one TRACOE phon assist I valve.

2. Product Description

TRACOE phon assist I valves were developed as adjustable speaking valves with or without an oxygen supply port for tracheostomy patients.

The head of the valve can be rotated on the base and contains a silicone membrane. The head and base sections each contain two lateral wedge-shaped apertures located opposite each other, which may be closed by rotating through 90 degrees clockwise and opened by turning anticlockwise. When the side apertures are closed, the valve operates as a one-way valve, i.e. it allows air to flow in only during inhalation, allowing voice production in the larynx during exhalation. After rotating to the left (anticlockwise), the patient can breathe in and out through the lateral valve apertures without obstruction. The lateral slot-shaped apertures are infinitely adjustable and can be opened wider or less wide as required.

In the base section of REF 650-SO and REF 650-TO/-TO-C, an oxygen supply port is provided in the form of a male Luer connector. This can either be connected to an oxygen supply using a connecting tube (REF 702) or closed using the occlusion cap supplied.

The phon assist I valves (REF 650-S and REF 650-SO) form an exact and air-tight seal with their base plates, both with the neck plate of the TRACOE grid and stoma buttons (REF 601 ff.) and with the cylindrical holder of the TRACOE adhesive carrier (REF 660).

The phon assist I valves (REF 650-T, REF 650-TO, REF 650-TO-C) each have a female adaptor for a 15 mm connector and therefore fit on all normal commercial tubes with 15 mm connectors.

3. Indications/Contraindications

3.1 Indications

TRACOE phon assist I valves are used for tracheotomized, awake patients with a completely or partly extant larynx and (at least partially present) deglutition abilities. They can be used in combination with a stoma button, grid button, adhesive carrier or a fenestrated tracheostomy tube with a 15 mm connector and allow the patient to speak. The valves can also be used when weaning the patient from a tracheostomy tube.

3.2 Contraindications

TRACOE phon assist I valves may under no circumstances be used in patients after complete removal of the larynx; this also applies to patients, who have been fitted with a voice prosthesis. Use of the valves when the tracheostomy tube is blocked is contraindicated.

If TRACOE phon assist I valves are used in combination with tracheostomy tubes, they should preferably be used with fenestrated tubes without a cuff. When used with fenestrated tubes with cuffs, the cuff must always be fully deflated.

CAUTION: TRACOE phon assist I valves must be removed when sleeping at night.

4. Use of TRACOE phon assist I Valves

4.1 Insertion and Removal of the Valve

Insertion of the REF 650-S and REF 650-SO into the TRACOE grid or stoma button should preferably be completed before introducing the button into the tracheostoma; in view of the simplicity and ease of this operation, the valve may, however, also be inserted in the

button when located in the tracheostoma. The special design of the valve base plate and the flexible silicone seal on the button guarantees a reliable and airtight fit between these parts. This also ensures that the valve is not detached from the button in the event of violent coughing.

REF 650-S and REF 650-SO may be inserted into and removed from the TRACOE adhesive carrier (REF 660) in a similar way.

REF 650-T, REF 650-TO and REF 650-TO-C are pushed onto the 15 mm connector of the fenestrated tracheostomy tube (e.g. TRACOE *twist*, TRACOE *twist plus* or TRACOE *comfort* tube).

The valve can be removed without difficulty from the button located in the tracheostoma or from the TRACOE *twist* cannula while still in place. The valve is removed either when the patient no longer wishes to speak or when secretions have collected in the valve and this needs to be cleaned.

4.2 General Requirements for Use of the Valves are as follows:

- adequate spontaneous breathing in the conscious patient
- stable cardiovascular conditions
- ability to swallow (at least partially present)
- no excessive bronchial secretion
- blood gas values in the normal range, with oxygen supply, if necessary

These requirements apply both when the valves are being used for speaking and when weaning the patient from a tracheostomy tube.

CAUTION: if, when speaking or during weaning, the patient becomes breathless or air becomes trapped, the lateral apertures are opened by turning the valve head anticlockwise as required. This disables the one-way valve function and the patient can breathe in and out without difficulty, without having to remove the valve.

4.3 Speaking with the Valve

For speaking, the lateral apertures of the valve are closed by rotating the valve head fully clockwise. This allows air to enter only during inspiration. During expiration, the valve closes and the air flows out through the larynx, throat and mouth, allowing voice production.

The continuous adjustability of the lateral apertures allows setting of an individual exhalation resistance. This can be used to increase oxygen saturation in the blood, if necessary.

4.4 Use of the Valves for Weaning

In patients, who have been ventilated with a tracheostomy tube for a longer period of time, weaning from the tube may be made easier by using a TRACOE phon assist I speaking valve. When weaning, the points mentioned in 4.2 must always be borne in mind.

Caution: the TRACOE phon assist I speaking valve must under no circumstances be used for weaning in combination with a blocked tracheostomy tube, as although the patient can still breathe in, breathing out is impossible. This would lead to hyperinflation of the lungs (air trapping) and possibly very soon to a life-threatening situation for the patient.

When using fenestrated tracheostomy tubes with a cuff, the cuff must be fully deflated before the speaking valve is fitted.

The enclosed warning label must be applied to the connector of the speaking valve and/or the inflation line of the tube to indicate the correct use of the speaking valve.

The possibility of opening the valve during weaning training means that air trapping can easily be avoided.

The patient must be monitored during use of the valve.

4.5 Further Information

If oxygen is to be supplied through valves REF 650-SO, REF 650-TO and REF 650-TO-C, care must always be taken to ensure that this is humidified and that the upper airways do not become dried out. The width of the lateral slots should also be taken into account here.

We recommend that at least one spare valve is kept available at all times.

If the function of the valve is impaired by mucus, the valve should be removed and cleaned.

5. Cleaning, Care and Disinfection

CAUTION: when cleaning the valves, care must be taken to ensure that the silicone membrane is not damaged.

5.1 Cleaning and Care

The speaking valve must be cleaned thoroughly after each use. Thoroughly rinse the valve with warm (but not hot) water and mild soap or detergent. Make sure that the speaking valve is fully dried. You may use a hair drier that is set to room temperature for drying.

CAUTION: TRACOE speaking valves must not be heated above 60 °C. The valve must not be cleaned with household or dental cleaning agents, hot water, peroxide, bleach, vinegar or alcohol, because these substances may impair the correct functioning of the valve. Do not heat up, boil or steam-sterilise the speaking valves.

5.2 Disinfection

Because the upper respiratory tract is never without germs – even in healthy individuals – we do not recommend the use of disinfectants.

If, however, disinfection appears necessary, TRACOE medical can provide, upon request, a list of tested disinfectants

The disinfected valve should then be rinsed very thoroughly with distilled water and dried.

CAUTION: Under no circumstances use disinfectants that release chlorine or substances that contain strong alkalis or phenol derivatives.

6. Duration of Use, Storage

Depending on how well it is looked after, the TRACOE phon assist I valve may be used for up to 6 months. We recommend that at least one spare valve is kept available. Even if valves are only slightly damaged, they must be replaced immediately. Unused valves should always be kept in a clean, dry condition in a clean ventilated container.

Storage: The products must be stored in a cool and dry location. Temporary storage for ≤ 30 days at $\leq +50$ °C is possible. Protect the product from sunlight and extreme cold.

CAUTION: before being used again, the valve must be examined to make sure that it is intact.

7. Returns

Returns of used products are accepted only if accompanied by a completed decontamination certificate. You can obtain the required form either directly from TRACOE medical or from the website www.tracoe.com.

8. General Terms and Conditions

The sale, delivery and return of all TRACOE products are carried out exclusively according to the General Terms and Conditions, which can be obtained directly from TRACOE medical GmbH or from the website www.tracoe.com.

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